# MCS2

# **USER MANUAL**









Copyright © 2018 SmarAct GmbH

Specifications are subject to change without notice. All rights reserved. Reproduction of images, tables or diagrams prohibited.

The information given in this document was carefully checked by our team and is constantly updated. Nevertheless, it is not possible to fully exclude the presence of errors. In order to always get the latest information, please contact our technical sales team.

SmarAct GmbH, Schuette-Lanz-Strasse 9, D-26135 Oldenburg Phone: +49 (0) 441 - 800879-0, Telefax: +49 (0) 441 - 800879-21

Internet: www.smaract.com, E-Mail: info@smaract.com

Document Version: 1.1.2

# TABLE OF CONTENTS

1	Maı	nufacturer Declarations	5
	1.1	Declaration of Conformity	5
	1.2	Information, Warnings and Safety Instructions	6
		1.2.1 Information – Intended Use	6
		1.2.2 Danger – Hazardous Voltage	6
		1.2.3 Caution – Unpacking	7
		1.2.4 Caution – Installation Instructions	7
		1.2.5 Caution – Connecting Instructions	8
		1.2.6 Caution – Disposal of Old Equipment	8
		1.2.7 Warranty and Liability	9
		1.2.8 Life Support Policy	9
2	Pro	duct Overview	10
	2.1	MCS2 Main Controller	11
		2.1.1 I/O Modules	13
	2.2	MCS2 Sensor Module	15
	2.3	MCS2 Accessories	15
3	Inst	tallation and Operation	16
	3.1	Setting up the Device	16
	3.2	Ethernet Interface Configuration	17
		3.2.1 Integrated Web Server	17
		3.2.2 Hand Control Module	19
		3.2.3 SmarActNetConfig Tool	19
	3.3	Finding An Unknown IP Address	19
4	Tecl	hnical Data	21
	4.1	Status LED	
	4.2	MCS2 Main Controller Connectors	
		4.2.1 Power Supply Connectors	
		4.2.2 D-SUB-15 Connector for Sensor Module	
		4.2.3 D-SUB-15-HD Connector for External Hand Control Module	
		4.2.4 USB Connector	
		4.2.5 RJ-45 Connector for Ethernet Communication	
		4.2.6 BNC Connectors for Fast Digital Input/Outputs	
		4.2.7 D-SUB-37 Connector of Auxiliary I/O	25
	4.3	MCS2 Sensor Module Connectors	
		4.3.1 D-SUB-15 Positioner Connectors	27
	4.4	MCS2 Electrical Specifications	28
		4.4.1 Positioner Driver Output	28
		4.4.2 Auxiliary Fast Digital Input/Outputs	28

## **TABLE OF CONTENTS**

4.4.3	Auxiliary GP Digital Inputs/Outputs	29
4.4.4	Auxiliary GP Analog Outputs	30
4.4.5	Auxiliary Fast Analog Inputs	31
4.4.6	Abbreviations	33

# **1** MANUFACTURER DECLARATIONS

# 1.1 Declaration of Conformity

according to DIN EN ISO/IEC 17050-1:2010

Manufacturer: SmarAct GmbH

Manufacturer's Address: Schuette-Lanz-Strasse 9

26135 Oldenburg, Germany



### The manufacturer hereby declares that the product

Product name: MCS2

Model Numbers: MCS2-C-xxxx, MCS2-S-xxxx, MCS2-H-xxxx,

MCS2-R19-xxxx, MCS2-MOD-xxxx

Product Options: al

complies – if installed in a compatible chassis from SmarAct and connected to a power supply provided by SmarAct – with the following European directives:

2014/35/EC Low Voltage Directive

2014/30/EC EMC Directive 2011/65/EU RoHS Directive

## The applied standards certifying the conformity are listed below:

Electromagnetic Emission: EN61000-6-3:2011, EN55011:2011 (Group 1, Class A)

Electromagnetic Immunity: **EN61000-6-1:2007**Safety (Low Voltage Directive): **EN61010-1:2011-07** 

July 3, 2018

Oldenburg, Germany

Managing Director

# 1.2 Information, Warnings and Safety Instructions

Please read the following information, warnings and safety instructions in this section carefully before using the product. The symbols and labels in this manual have the following definitions:

# Symbol Definition



This symbol identifies conditions or practices that could result in damage to the equipment or other property, and in extreme cases, possible personal injury or death.



This symbol identifies conditions or practices that involve potential electric shock hazard.

#### 1.2.1 Information - Intended Use

The SmarAct MCS2 controllers are designed to drive SmarAct positioners only. Never connect any other equipment than SmarAct positioners to the MCS2 controller outputs.

SmarAct MCS2 controllers are laboratory equipment designed according to the safety requirements for electrical equipment for measurement, control and laboratory use (DIN EN 61010).

The MCS2 controller should be operated under the following environmental conditions:

- Indoor use only
- Altitude up to 2000 m
- Temperature range from 5 °C to 40 °C
- Maximum relative humidity 80% up to 31  $^{\circ}$ C, decreasing linearly to 50% relative humidity at 40  $^{\circ}$ C
- Degree of pollution: 2

Please use the provided power supply only. This is required to meet all design specifications.

# 1.2.2 Danger - Hazardous Voltage

The MCS2 controller described in this manual is capable of generating high output currents at high voltages. They may cause serious or even lethal injury if used improperly. Therefore, the equipment should only be operated by personnel that is adequately trained and educated to prevent any improper use.

Please note the following general accident prevention rules:



### **DANGER**

- Never touch any part that might be connected to an output with a high voltage.
- Do not connect products from other manufacturers to the output connectors.
- · Never use equipment that is damaged in any way.
- The MCS2 controller contains no user serviceable parts. Never open the housing. Procedures which require to open the housing must only be carried out by authorized, qualified and trained personnel.

Output connectors with dangerous signals are labeled with the following symbol:





#### WARNING

The Sensor Modules connected to these output connectors have pins with a high voltage, too.

# 1.2.3 Caution - Unpacking

Please be careful when unpacking the MCS2 controller. Inspect the MCS2 controller for signs of damage and only use equipment that shows no signs of damage. In case of any damage, contact SmarAct for replacement. Please save all packing materials in case you would like to transport or ship the product again.

#### 1.2.4 Caution - Installation Instructions

The MCS2 controller must be installed horizontally with 3 cm air circulation area behind the fan. Insufficient air flow can cause overheating, which can result in a limited functionality of the controller.

If you are using an OEM version of an MCS2 controller, the controller board must be installed in such a way, that free air convection can be achieved. Insufficient ventilation will cause overheating and premature failure. Therefore forced air cooling, for example in a rack with a fan, is preferred. Make sure that no part of the MCS2 OEM controller board other than the required connectors can be touched after final installation. Please take care to meet all applicable requirements concerning electrical equipment in your regulatory domain.



### **CAUTION**

Avoid touching any circuit components, pins or PCB traces, as these components are ESD-sensitive devices.

Please touch the OEM boards on their edges only during handling and installation. Discharge any static charge on your body by touching a grounded object prior to the handling of the MCS2 controller board. If not installed, place the MCS2 controller boards only on conductive surfaces, such as ESD- safe transport containers (envelopes or foam).

Only apply power to the MCS2 OEM controller board if it is installed in the intended housing. Never apply power to a MCS2 OEM controller board if any other part than the required connectors can be touched.

### 1.2.5 Caution - Connecting Instructions

Never use any other connecting cables than the connecting gear that you obtained from SmarAct. Never use any third party adapters or cables. This can cause failure or malfunction.



#### NOTICE

When removing the USB or Ethernet cable, all positioners will be stopped immediately as a safety precaution.

# 1.2.6 Caution – Disposal of Old Equipment

According to the EU directive 2012/19/EU, as of August 13th, 2012, electrical and electronic equipment may not be disposed of with other wastes in the member states of the European Union.

SmarAct respects the manufacturer's product responsibility and will take care of environmentally correct disposal of old SmarAct products free of charge.

To dispose of old SmarAct equipment, you can return it to SmarAct to the following address postage-free:

SmarAct GmbH Schuette-Lanz-Strasse 9 26135 Oldenburg, Germany



WEEE-Reg.-Nr.: DE 47992153

### 1.2.7 Warranty and Liability

The General Terms and Conditions of Sale and Delivery from SmarAct GmbH always apply. These conditions are available to the operator upon signing the contract, placing an order or at http://www.smaract.com/site-notice/. No warranty or liability claims may be made in the event of injury to persons or damage to property if this has arisen from one or more of the following:

- · Improper use of the unit,
- · incorrect assembly, startup and operation,
- · operation with defective and/or non-functional safety and protective devices,
- failure to comply with the information in the operating instructions regarding safety, transport, storage, mounting/installation, commissioning, operation and maintenance of the unit,
- · incorrectly or unauthorized repairs,
- unauthorized modifications to the device.
- inadequate monitoring of parts which are subject to wear,
- damage caused by exposure to water, e.g. condensation water formation, where this is outside the responsibility of SmarAct GmbH,
- · the effect of foreign bodies or mechanical damage,
- · disastrous occurrences,
- · force majeure.

Please refer to section 1.2.6 for the contact address for claims under warranty, repair and replacement service.

### 1.2.8 Life Support Policy

SmarAct GmbH does not authorize or warrant any of its products for use in life support systems, without the specific written consent of SmarAct GmbH.

Life support systems are equipment intended to support or sustain life, and whose failure to perform, when properly used in accordance with instructions provided, can be reasonably expected to result in personal injury or death.

# **2** PRODUCT OVERVIEW

This document is a user manual for the SmarAct Modular Control System 2 (MCS2) which is designed to drive piezo based SmarAct stepping actuators.

For the general concepts of controlling the MCS2 by software as well as for a description of the application programming interface (API) please refer to the MCS2 Programmers Guide.

Each MCS2 control device consists of the following components:

### MCS2 Main Controller

The main controller has the following tasks:

- process commands from the hand control module or from the user application software
- read in the sensor data from the integrated position sensors (via sensor module)
- perform closed-loop position control
- drive the positioners

Each MCS2 controller is equipped with a USB or Network Interface and can be controlled by software running on a PC.

### MCS2 Sensor Module

The positioners may be equipped with integrated sensors to perform closed-loop positioning control. For this, a sensor module preprocesses the sensor data and delivers the digitized data to the main controller. Alternatively, a SmarAct PicoScale interferometer can be used as a sensor module. Each sensor module can handle up to three positioners.

### Hand Control Module (optional)

The MCS2 offers easy and ready-to-go control by joysticks and control knobs - without the need of complex installation procedures or software development. The Hand Control Module may be integrated into the main controller or placed inside a separate (external) housing.

#### I/O Module (optional)

The MCS2 can be equipped with I/O modules that provide digital and analog inputs and outputs. Fast digital outputs can be used to trigger external devices on specific internal events of the main controller (e.g. Position Reached). Digital device inputs allow to trigger the MCS2 by an external device for example to perform an Emergency Stop or to synchronize data streaming. General purpose digital inputs and outputs provide control signals to switch lights, relays, dispensers, etc. or to read the state of safety switches, light barriers, etc. Analog inputs can be used to read analog voltage signals from external devices. The controller supports to feed these signals into the control loop. Closed loop operation depending on external sensor signals can be implemented with this feature.

The I/O module is integrated into the main controller.

Figure 2.1 shows a structural setup of a typical MCS2 system.




Figure 2.1: Structural Setup of an MCS2 Device

Note that all external connections support hot-plugging. The controller will automatically detect and report a changed configuration. If the connection to the PC or Hand Control Module is unplugged, all movements that were sent over the disconnected interface are stopped as a safety measure.

As there is a wide range of applications, many different combinations of MCS2 main controllers and sensor modules are configurable to meet the individual needs of our customers. Sections 2.1 and 2.2 give an overview on current models.

### 2.1 MCS2 Main Controller

MCS2 main controllers are available with different housing, a selectable number of channels and an interface option which offers all the connectivity needed. The configuration of your product is defined by the product code below and the options listed in table 2.1.




Table 2.1 - Product Codes



Continued on next page

Table 2.1 – Continued from previous page

Product Code	Description	Connectors
MCS2-C-0001	3 channel Controller, integrated Hand Control Module, USB interface	USB type B, D-SUB-15 female (Sensor Module)
MCS2-C-0002	3 channel Controller, USB interface	USB type B, D-SUB-15 female (Sensor Module), D-SUB-15-HD female (ext. Hand Control Module)
MCS2-C-0003	6 channel Controller, USB interface	USB type B, D-SUB-15 female (Sensor Module)
MCS2-C-0004	9 channel Controller, USB interface	USB type B, D-SUB-15 female (Sensor Module)
MCS2-C-0005	3 channel Controller, integrated Hand Control Module, Ethernet interface	RJ-45 Ethernet, D-SUB-15 female (Sensor Module)
MCS2-C-0006	3 channel Controller, Ethernet interface	RJ-45 Ethernet, D-SUB-15 female (Sensor Module)
MCS2-C-0007	6 channel Controller, Ethernet interface	RJ-45 Ethernet, D-SUB-15 female (Sensor Module)
MCS2-C-0008	9 channel Controller, Ethernet interface	RJ-45 Ethernet, D-SUB-15 female (Sensor Module)
MCS2-C-0009	3 channel Controller, I/O Module Type 3, USB interface	USB type B, D-SUB-15 female (Sensor Module), 4x BNC (3 out, 1 in)
MCS2-C-0010	12 channel Controller, USB interface	USB type B, D-SUB-15 female (Sensor Module)
MCS2-C-0011	12 channel Controller, Ethernet interface	RJ-45 Ethernet, D-SUB-15 female (Sensor Module)
MCS2-C-0012	3 channel Controller, I/O Module Type 5, USB interface	USB type B, D-SUB-15 female (Sensor Module), 3x BNC (3 out), USB type B (HSDR)
MCS2-C-0013	3 channel Controller, I/O Module Type 5, Ethernet interface	RJ-45 Ethernet, D-SUB-15 female (Sensor Module), 3x BNC (3 out), USB type B (HSDR)
MCS2-C-0014	6 channel Controller, I/O Module Type 3, Ethernet interface	RJ-45 Ethernet, D-SUB-15 female (Sensor Module), 4x BNC (3 out, 1 in)
MCS2-C-0015	6 channel Controller, I/O Module Type 1/3, Ethernet interface	RJ-45 Ethernet, D-SUB-15 female (Sensor Module), 7x BNC (6 out, 1 in)
MCS2-C-0016	3 channel Controller, I/O Module Type 3, Ethernet interface	RJ-45 Ethernet, D-SUB-15 female (Sensor Module), 4x BNC (3 out, 1 in)

Continued on next page

Table 2.	l – Continuea	from previous <sub>l</sub>	page
_			

Product Code	Description	Connectors
MCS2-C-0017	3 channel Controller, I/O Module Type 4, USB interface	USB type B, D-SUB-15 female (Sensor Module), D-SUB-37 female (IO) 1x BNC (1 in)
MCS2-C-0018	3 channel Controller, I/O Module Type 4, Ethernet interface	RJ-45 Ethernet, D-SUB-15 female (Sensor Module), D-SUB-37 female (IO) 1x BNC (1 in)
MCS2-C-0019	9 channel Controller, I/O Module Type 3, USB interface	USB type B, D-SUB-15 female (Sensor Module), 4x BNC (3 out, 1 in)
MCS2-C-0020	9 channel Controller, I/O Module Type 3, Ethernet interface	RJ-45 Ethernet, D-SUB-15 female (Sensor Module), 4x BNC (3 out, 1 in)

The connectors of the most common MCS2 controllers are shown in the following figures.




Figure 2.2: Three Channel MCS2 Table Top Housing (front view)

Depending on the computer interface the arrangement is a bit different. Controllers with more channels that are delivered in a table-top housing have the same connectors as in the figures, but are equipped with additional D-SUB-15 connectors to the sensor modules.

### 2.1.1 I/O Modules

Depending on a specific product code the MCS2 may be equipped with one or more I/O modules. The following tables list the currently available I/O module types and their features:

### I/O Module Type 1

Table 2.2 – I/O Module Type 1

Name	Туре	Connectors	Function
DOUT-n	3 Fast Digital Outputs	BNC	Output Trigger: Position Compare, Position Reached, Actively Moving

# I/O Module Type 3

Table 2.3 – I/O Module Type 3

Name	Туре	Connectors	Function
DIN-1	1 Digital Device Input	BNC	Input Trigger: Emergency Stop, Trajectory Streaming Synchronization, Command Group Trigger
DOUT-n	3 Fast Digital Outputs	BNC	Output Trigger: Position Compare, Position Reached, Actively Moving

# I/O Module Type 2

Table 2.4 – I/O Module Type 2

Name	Туре	Connectors	Function
DOUT-n	3 Fast Digital Output	D-SUB37	Output Trigger: Position Compare, Position Reached, Actively Moving
GP-DIN-n	4 Digital Inputs	D-SUB37	General Purpose Digital Input
GP-DOUT-n	4 Digital Outputs	D-SUB37	General Purpose Digital Output, Open- Collector
AOUT-n	2 Analog Outputs	D-SUB37	General Purpose +/-10V Analog Voltage Output
AIN-n	6 Fast Analog Inputs	D-SUB37	General Purpose Analog Input, Control Loop Feedback

# I/O Module Type 4

Table 2.5 – I/O Module Type 4

Name	Туре	Connectors	Function
DIN-1	1 Digital Device Input	D-SUB37 / BNC	Input Trigger: Emergency Stop, Trajectory Streaming Synchronization, Command Group Trigger

Continued on next page

table 2.3 Continued from previous page					
Name	Туре	Connectors	Function		
DOUT-n	3 Fast Digital Output	D-SUB37	Output Trigger: Position Compare, Position Reached, Actively Moving		
GP-DIN-n	4 Digital Inputs	D-SUB37	General Purpose Digital Input		
GP-DOUT-n	4 Digital Outputs	D-SUB37	General Purpose Digital Output, Open- Collector		
AOUT-n	2 Analog Outputs	D-SUB37	General Purpose +/-10V Analog Voltage Output		
AIN-n	6 Fast Analog Inputs	D-SUB37	General Purpose Analog Input, Control		

Table 2.5 - Continued from previous page

## 2.2 MCS2 Sensor Module

For positioners with integrated sensors a sensor module is required to digitize the sensor data. The sensor module type depends on the sensor type, the driver type and the connector to the positioner. The configuration of the sensor module is defined by the product code below. In table 2.6 you will find all available configuration options.




Table 2.6 - Product Codes

Product Code	Description	Connectors
MCS2-S-0001	Sensor Module for up to 3 positioners	3x D-SUB-15 female

## 2.3 MCS2 Accessories

In addition to the modules above, some optional accessories may be required for specific configurations. In table 2.7 you will find a list of available accessories.

Table 2.7 - Product Codes

Product Code	Description	Connectors
MCS2-A-PS-CABLE-1.5M-1.5M	Adapter cable for MCS2 and PicoScale	3x D-SUB-15 female

# **3** INSTALLATION AND OPERATION

# 3.1 Setting up the Device

Before initial operation be sure to read the safety notes in chapter 1. To set up your MCS2 device please follow the instructions below:

- 1. Make sure that the **power switch** of your MCS2 controller is in the **off** position (0) or that the MCS2 control module is not supplied with power, respectively.
- 2. **Connect** the **Sensor Modules**, which may be integrated into the manipulator base plate, to the D-SUB 15 pin female connector(s) of the MCS2 controller. If the connectors are labeled, their labels must match.
- 3. If the cables from the positioners or Sensor Modules are delivered with feed-throughs, **install the feed-throughs** at a flange of the vacuum chamber and connect the matching connectors to both sides of the feed-through. If the connectors are labeled, their labels must match.
- 4. If the Sensor Modules are delivered in a separate housing, **connect the positioners** to the D-SUB 15 pin female connectors of the Sensor Modules. If the connectors are labeled, their labels must match.
- 5. If the Hand Control Module is delivered in a separate housing, connect the **Hand Control Module** to the D-SUB 15 pin high-density female connector of the MCS2 main controller using the enclosed cable.
- 6. If desired, **connect the ground potential** of the mechanical system setup to the ground pin of the MCS2 controller. This is recommended for applications in scanning electron microscopes where the ground potential of the system setup should be equivalent to the ground potential of the controllers.
- 7. **Connect the power supply**, which has been delivered with the controller, to the power supply connector of the MCS2 controller.
- 8. **Connect the MCS2 controller** with one of the following **interface** cables:
  - · USB cable to the PC
  - Ethernet cable to a switch or PC
- 9. Power up the controller by switching the **power switch** to the **on** position (I).
- 10. The positioning system now can be controlled by software. Make sure that the correct **positioner types** are configured for your positioners and perform the **calibration routine** on the controller for each channel.

For getting started the *MCS2 Demo GUI* delivered with each MCS2 device may be used. The *MCS2 Service Tool* allows to easily configure and setup the MCS2 controller. For detailed information on these programs please refer to the specific manuals.

To create your personal software SmarAct delivers flexible and well documented Windows® and Linux libraries. For those who like to use LabVIEW™, we provide ready-to-use VIs based on the Windows® DLL. Programming examples for C and LabVIEW™allow the developer a quick start into the application programming with the MCS2.

For more information on the available software please refer to the MCS2 Software Installation documentation.

# 3.2 Ethernet Interface Configuration

While devices with USB interface do not need any interface configuration, the ethernet interface must be configured with the network parameters: DHCP mode, IP address, subnet mask and gateway IP address. The MCS2 is delivered with a default IP configuration which may be adjusted to match the users network settings.

The following table lists the default configuration:

Parameter	Default Value		
DHCP Mode	disabled		
IP Address	192.168.1.200		
Subnet Mask	255.255.0.0		
Gateway IP	192.168.1.1		
Pass-Key	smaract		

The MCS2 offers several ways to configure the network parameters. These are described in detail in the following sections.

In general, the interface may be configured to use DHCP to obtain an IP address from a DHCP-server or to use a static IP configuration.

### 3.2.1 Integrated Web Server

The built-in web server is accessible with any internet browser (e.g. Firefox®). Obviously, the current IP address of the controller must be known to be able to connect to the web server. Furthermore, the MCS2 and the PC need to be on the same network to enable the browser to load the configuration web page.

- Open the internet browser and type the IP address of the controller in the address bar. The default IP is: 192.168.1.200.
- Press 'Enter' to access the webserver.

The following webpage should be displayed:




Figure 3.1: Device Information Tab of the Webpage

The 'Information' tab gives some general information about the device, like serial number, device name, number of channels, etc.




Figure 3.2: Configuration Tab of the Webpage

The **'Configuration'** tab shows the current network configuration and allows to configure the interface. Note that changes only become effective after the configuration was saved. The configuration is protected against unintentional modification by a pass-key.

- To enable the **DHCP** mode, select the 'DHCP' option. No further settings are required in that case. After saving the configuration, the controller will try to get an IP address automatically from a DHCP server. If no server is available for a certain period of time, the interface will fall-back to the static default IP settings. Click the 'Save Configuration' button, enter the pass-key ('smaract') in the upcoming dialog window and click 'Ok'.
- To configure a **static IP**, select 'Static IP' and enter the IP address, the gateway IP and the subnet mask in the designated input fields. Click the 'Save Configuration' button, enter the pass-key ('smaract') and click 'Ok'.

• To **reset** the interface to the factory **default** settings click the 'Reset to Default' button, enter the pass-key ('smaract') and click 'Ok'.



#### NOTICE

The configuration is saved permanently and the interface will restart with the new configuration! Make sure that the given parameters match your network setup, otherwise the configuration page may become unreachable!

#### 3.2.2 Hand Control Module

The Hand Control Module of the MCS2 may be used to configure the network interface settings. Enter the 'Main Menu' -> 'Settings' -> 'Network Config.' menu and modify the parameters as desired. Select 'Save Config.' to apply the changes.

### 3.2.3 SmarActNetConfig Tool

The **SmarActNetConfig** tool may be used to find and configure MCS2 devices.

The tool is installed with the MCS2 Software Installer. See the *SmarAct NetConfig User Manual* document for more details on the tool.

# 3.3 Finding An Unknown IP Address

There are two methods to find an unknown IP address of a MCS2 controller:

- 1. Use the **SmarActNetConfig** tool
- 2. Set the MCS2 to default configuration with the 'default config' **push button**

The **SmarActNetConfig** tool is able to find devices regardless of its actual network configuration, even if the parameters do not match the current network settings.

Note that the device cannot be found if the discovering was disabled with the *Network Discover Mode* property.

For the second method the MCS2 device has a **push button** on the front panel which may be used to reset the interface configuration to the static default configuration temporarily. The button is located behind a small hole below the ethernet connector on the front panel. Use a pen or a paper-clip to press the button for minimum 5 seconds. Note that this IP setting is volatile, meaning that a power cycle will restore the original settings assumed they have not been reconfigured in this session.

Now connect to the MCS2 configuration webpage with an internet browser by using the default IP address: 192.168.1.200. Use the configuration tab to configure a new IP address as desired. For this method the MCS2 must be connected to a network in the address range: 192.168.xxx.xxx. The

easiest way to do this is to connect the MCS2 directly to the PC ethernet interface and to configure a static IP address for the PC. (For example 192.168.1.100).

# 4 TECHNICAL DATA

### 4.1 Status LED

MCS2 controllers without an integrated Hand Control Module are equipped with a status LED. This LED shows the current device state. Usually the LED is permanently on signalizing the normal operation. A sensor module attach or detach process is indicated by fading the LED down and up for ca. 1 second. According to the following table different flashing patterns may occur.

LED Pattern	Frequency	Description
On	-	Regular operation
Flashing	2 Hz	Over-temperature of min. one amplifier
Flashing	5 Hz	High-voltage overload
Flashing	10 Hz	High-voltage supply failure
Flashing	15 Hz	The device is in the bootloader state and has no valid application
Fading up and down	1 Hz	The device is in the bootloader state

Table 4.1 – LED states

# 4.2 MCS2 Main Controller Connectors

## **4.2.1 Power Supply Connectors**

The MCS2 controller is supplied by a 12V power supply. Please use the power supply that has been shipped with the MCS2 controller and make sure it matches the current and power rating, which is specified on the label of the controller (on the bottom or back side of the housing).

#### 3-/6-Channel Controller



Figure 4.1: DC power connector for power supply

Table 4.2 – DC Power Connector Pin Assignment

Pin	Signal	Function		
Tip	12V	+12V DC Controller Power Supply		
Ring	RTN	Power Supply Return		

## 9-Channel Controller




Figure 4.2: Female KPPX connector for power supply

Table 4.3 - KPPX Pin Assignment

Pin	Signal	Function		
1, 2	12V	+12V DC Controller Power Supply		
3, 4	RTN	Power Supply Return		
Shielding	PE	Protective Earth		

### **12-Channel Controller**





Figure 4.3: Female Mini Fit Jr. connector for power supply

Table 4.4 – Mini Fit Jr. Pin Assignment

Pin	Signal	Function		
1, 2, 3	12V	+12V DC Controller Power Supply		
4, 5, 6	RTN	Power Supply Return		

## 4.2.2 D-SUB-15 Connector for Sensor Module

The MCS2 controller provides one or more female D-SUB-15 connectors to stick-slip positioners via an MCS2 Sensor Module. The pin layout is described below.




D-SUB-15 connector to Sensor Module

Table 4.5 – DSUB15 Pin Assignment

Pin	Signal	Function
1	HV-OUT-1	Positioner driving signal, channel 1, 4, 7,
2	HV-OUT-2	Positioner driving signal, channel 2, 5, 8,
3	HV-OUT-3	Positioner driving signal, channel 3, 6, 9,
4	SM-GND	Ground for Sensor Module and sensor
5	SM-TX-D-	RS-422 D- signal from Sensor Module
6	SM-RX-D-	RS-422 D- signal to Sensor Module
7	d.n.c.	DO NOT CONNECT
8	EXT_STATE	Internal signal lines
9	HV-GND-1	Ground for positioner driving signal, channel 1, 4, 7,
10	HV-GND-2	Ground for positioner driving signal, channel 2, 5, 8,
11	HV-GND-3	Ground for positioner driving signal, channel 3, 6, 9,
12	SM-TX-D+	RS-422 D+ signal from Sensor Module
13	SM-RX-D+	RS-422 D+ signal to Sensor Module
14	SyncB	Internal sync line
15	SM-5V	Power supply for Sensor Module, 5V DC
Shielding	SM-Shield	Shielding for Sensor Module and sensor

## 4.2.3 D-SUB-15-HD Connector for External Hand Control Module

For external Hand Control Modules the MCS2 controller provides a female D-SUB-15 high-density connector. The pin layout is described below.



D-SUB-15-HD connector to external Hand Control Module

Table 4.6 – Pin Assignment

Pin	Signal	Function
1	GND	Ground
2	d.n.c.	DO NOT CONNECT
3	d.n.c.	DO NOT CONNECT
4	d.n.c.	DO NOT CONNECT
5	d.n.c.	DO NOT CONNECT
6	3.3V	Power supply, 3.3V DC
7	d.n.c.	DO NOT CONNECT
8	GND	Ground
9	SyncB	Internal sync line
10	5V	Power supply, 5V DC
11	Switch	Power switch
12	RXD+	Differential data for receiving
13	RXD-	Differential data for receiving
14	TXD+	Differential data for transmitting
15	TXD-	Differential data for transmitting
Shielding	GND	Ground

#### 4.2.4 USB Connector


The USB Type B connector offers direct connectivity e.g. to a PC for software control of the MCS2. The USB interface is galvanically isolated to the internal ground of the controller.

# 4.2.5 RJ-45 Connector for Ethernet Communication

The RJ-45 connector at the MCS2 main controller provides a standard 10BaseT / 100BaseTX ethernet interface for software control of the MCS2.

## 4.2.6 BNC Connectors for Fast Digital Input/Outputs


The BNC connectors provide fast digital outputs (DO1, DO2, DO3) and digital device inputs (DI1). The connectors are galvanically isolated to the internal ground of the controller.



**BNC** connector

# 4.2.7 D-SUB-37 Connector of Auxiliary I/O

The D-Sub-37 female connector provides analog and digital auxiliary inputs and outputs. All signals are galvanically isolated to the internal ground of the controller.



D-SUB-37 connector of Auxiliary I/O

Table 4.7 – Pin Assignment

Pin	Signal	Function
25	DIN-1	Digital Device Trigger Input
24	GNDiso	Corresponding GNDiso for DIN-1
21	GP-DIN-1	General Purpose Digital Input 1
22	GP-DIN-2	General Purpose Digital Input 2
23	GP-DIN-3	General Purpose Digital Input 3
6	GP-DIN-4	General Purpose Digital Input 4
20	GNDiso	Corresponding GNDiso for GP-DIN-1/-2/-3/-4
7	GP-DOUT-1	General Purpose Digital Output 1, Open Collector
8	GP-DOUT-2	General Purpose Digital Output 2, Open Collector
27	GP-DOUT-3	General Purpose Digital Output 3, Open Collector
9	GP-DOUT-4	General Purpose Digital Output 4, Open Collector
26	GNDiso	Corresponding GNDiso for GP-DOUT-1/-2/-3/-4
10	DOUT-1	Fast Digital Output 1
29	DOUT-2	Fast Digital Output 2
11	DOUT-3	Fast Digital Output 3
28	GNDiso	Corresponding GNDiso for DOUT-1/-2/-3


Continued on next page

Table 4.7 – Continued from previous page

Pin	Signal	Function
12	AIN-6	Fast Analog Input 6
30	GNDiso	Corresponding GNDiso for AIN-6
13	AIN-5	Fast Analog Input 5
31	GNDiso	Corresponding GNDiso for AIN-5
14	AIN-4	Fast Analog Input 4
32	GNDiso	Corresponding GNDiso for AIN-4
15	AIN-3	Fast Analog Input 3
33	GNDiso	Corresponding GNDiso for AIN-3
16	AIN-2	Fast Analog Input 2
34	GNDiso	Corresponding GNDiso for AIN-2
17	AIN-1	Fast Analog Input 1
35	GNDiso	Corresponding GNDiso for AIN-1
18	AOUT-1	General Purpose Analog Output 1
36	GNDiso	Corresponding GNDiso for AOUT-1
19	AOUT-2	General Purpose Analog Output 2
37	GNDiso	Corresponding GNDiso for AOUT-2
1	d.n.c	DO NOT CONNECT, Internally Used
2-5	NC	Not Connected
Shielding	GND	Ground

# **4.3 MCS2 Sensor Module Connectors**

## **4.3.1 D-SUB-15 Positioner Connectors**



Sensor Module D-SUB-15 connectors

Table 4.8 – Pin Assignment

Pin	Signal	Function
1	HV-OUT-x	Positioner driving signal, channel x
2	d.n.c.	DO NOT CONNECT
3	d.n.c.	DO NOT CONNECT
4	S-GND	Ground for sensor
5	S-SIN+	sin+ signal from sensor
6	S-COS+	cos+ signal from sensor
7	S-REF+	ref+ signal from sensor
8	d.n.c.	DO NOT CONNECT
9	HV-GND-x	Ground for positioner driving signal, channel x
10	d.n.c.	DO NOT CONNECT
11	d.n.c.	DO NOT CONNECT
12	S-SIN-	sin- signal from sensor
13	S-COS-	cos- signal from sensor
14	S-REF-	ref- signal from sensor
15	S-VCC	Power supply for sensor

\*The HV-OUT-x signals are identical to the HV-OUT-x signals from the MCS2 controller.

# 4.4 MCS2 Electrical Specifications

# **4.4.1 Positioner Driver Output**

Table 4.9 - Driver Output Specification

Parameter	Value	Unit
Driving Signal (HV-OUT-x)		
Signalform	sawtooth (step mode) constant (scan mode)	
Output voltage range	0 to 100	V
Average current per channel	200	mA
Peak current per channel, $<$ 10 $\mu$ s, max. speed	20	Α

# 4.4.2 Auxiliary Fast Digital Input/Outputs

All signals are galvanically isolated to the internal ground of the controller. The **digital device input** has schmitt-trigger characteristics. Table 4.10 lists the threshold voltage levels. The output voltage level of the **fast digital outputs** can be configured to 3.3V or 5V by software. This setting is global for all digital outputs. Note that the digital outputs are in a high-impedance state by default, therefore the output driver must be explicitly enabled.

Table 4.10 – Fast Digital Input/Output Specification

Parameter	Condition	Vcc	MIN	TYP	MAX	Unit		
	Common							
$V_{iso}$					60	V		
R <sub>iso</sub>				10		МΩ		
	Dig	ital Device i	Inputs					
$V_i$			-0.3		5.5	V		
$V_{ih}$		3.3V/5V		1.9	2.1	V		
$V_{il}$		3.3V/5V	1.0	1.2		V		
$R_i$				1		МΩ		
Fast Digital Outputs								
Io		3.3V			8	mA		
lo		5V			3	mA		

Continued on next page

Table 4.10 – Continued from previous page

Parameter	Condition	Vcc	MIN	TYP	MAX	Unit
V <sub>oh</sub>	Ioh = 20uA	3.3V/5V	Vcc-0.1			V
$V_{oh}$	Ioh = 2mA	3.3V	3.05			V
$V_{oh}$	Ioh = 3mA	3.3V	2.90			V
$V_{oh}$	Ioh = 5.5mA	3.3V	2.60			V
$V_{oh}$	Ioh = 8mA	3.3V	2.50			V
$V_{oh}$	Ioh = 2mA	5V	4.75			V
$V_{oh}$	Ioh = 3mA	5V	4.60			V
$V_{oh}$	Ioh = 5.5mA	5V	4.30			V
$V_{oh}$	Ioh = 8mA	5V	4.20			V
V <sub>ol</sub>	Iol = -20uA	3.3V/5V			0.10	V
V <sub>ol</sub>	Iol = -2mA	3.3V			0.20	V
V <sub>ol</sub>	Iol = -3mA	3.3V			0.26	V
V <sub>ol</sub>	Iol = -5.5mA	3.3V			0.49	V
V <sub>ol</sub>	Iol = -8mA	3.3V			0.71	V
V <sub>ol</sub>	Iol = -2mA	5V			0.20	V
V <sub>ol</sub>	Iol = -3mA	5V			0.25	V
V <sub>ol</sub>	Iol = -5.5mA	5V			0.48	V
V <sub>ol</sub>	Iol = -8mA	5V			0.70	V

# **4.4.3 Auxiliary GP Digital Inputs/Outputs**

All signals are galvanically isolated to the internal ground of the controller. The **general purpose digital inputs** have schmitt-trigger characteristics. Table 4.11 lists the threshold voltage levels. The **general purpose digital outputs** are designed as open-collector outputs. This means that the output logic is inverted. Writing a one to an output switches the output transistor on which leads to a low signal level at the output pin.

Table 4.11 – General Purpose Digital Inputs/Outputs Specification

Parameter	Condition	Vcc	MIN	TYP	MAX	Unit
	(	Common				
$V_{iso}$					60	V
R <sub>iso</sub>				10		МΩ
General Purpose Digital Inputs						

Continued on next page

Parameter	Condition	Vcc	MIN	TYP	MAX	Unit
$V_i$			-0.3		5.5	V
$V_{ih}$		3.3V/5V		1.9	2.1	V
$V_{il}$		3.3V/5V	1.0	1.2		V
R <sub>i</sub>				1		МΩ
	General Purp	pose Digita	l Output	ts		
l <sub>ol</sub>	Open Collector				700	mA
l <sub>oh</sub>	Open Collector			2		μΑ
V <sub>oh</sub>	Open Collector				24	V

Table 4.11 – Continued from previous page

# 4.4.4 Auxiliary GP Analog Outputs

 $R_{on}$ 

The **analog outputs** feature a DAC with 16bit resolution. All outputs are galvanically isolated to the internal ground of the controller. The corresponding GND signals of the specific outputs should be used to archive the best performance.

Open Collector



### NOTICE

Note that the analog outputs may generate a short random voltage peak in the range of  $\pm 10$ V while powering up or down the controller before the level stabilizes at its default level of 0V. External equipment must be capable of standing this or must be disconnected while power-up / power-down.

400

 $m\Omega$ 

Table 4.12 – General Purpose Analog Outputs Specification

Parameter	Condition	MIN	TYP	MAX	Unit		
Common							
$V_{iso}$				60	V		
R <sub>iso</sub>			10		ΜΩ		
	General Purpose Analog Outputs						
$V_0$		-10		10	V		
Io				10	mA		
R <sub>o</sub>			10		Ω		
INL	Using line passing through codes 512 and 65024		±4		LSB		
DNL			±0.2		LSB		

Continued on next page

Table 4.12 – Continued from previous page

Parameter	Condition	MIN	TYP	MAX	Unit
$V_{\text{offset}}$			±4		mV
$V_{\text{noise}}$			860		$\mu V_{pp}$
t <sub>settling</sub>			290		μs
Slew rate	measured btw. 10 to 90 % of a full-scale transition		0.163		V/µs
DC crosstalk	channel to channel		60		μV

# **4.4.5 Auxiliary Fast Analog Inputs**

The **fast analog inputs** feature an ADC with 16bit resolution and a build-in 2nd-order low pass filter. All inputs are galvanically isolated to the internal ground of the controller. The corresponding GND signals of the specific inputs should be used to archive the best performance. The input range of the analog inputs is configurable for unipolar as well as bipolar operation. Note that this setting is global for all analog inputs. The total number of six analog inputs are mapped in groups of two to the channels of the corresponding driver module.

Table 4.13 – Fast Analog Inputs Specification

Parameter	Condition	MIN	TYP	MAX	Unit
	Commo	on			
$V_{iso}$				60	V
R <sub>iso</sub>			10		МΩ
	Fast Analog	Inputs			
V <sub>ovp</sub>		-11		11	V
$R_i$			500		kΩ
$F_{co}$	-3dB		15		kHz
F <sub>sample</sub>	per channel		50		kHz
INL			±1	2	LSB
DNL			±0.5	1.5	LSB
SINAD			91.5		dB
Egain			0.1		% FSR
E <sub>offset</sub>			±0.5		mV
$V_i$	input range=5.12V	0		5.12	V
Vi	input range=10.24V	0		10.24	V
$V_i$	input range=±2.56V	-2.56		2.56	V

Continued on next page

Table 4.13 – Continued from previous page

Parameter	Condition	MIN	TYP	MAX	Unit
$V_i$	input range=±5.12V	-5.12		5.12	V
$V_i$	input range=±10.24V	-10.24		10.24	V

#### 4.4.6 Abbreviations

- V<sub>iso</sub>: allowed voltage difference from isolated ground to internal ground
- R<sub>iso</sub>: isolation resistance from isolated ground to internal ground
- V<sub>cc</sub>: configured common supply voltage
- V<sub>oh</sub>: minimal output value if output is set high
- Vol: maximal output value if output is set low
- V<sub>o</sub>: output voltage range
- I<sub>o</sub>: allowed output current
- I<sub>oh</sub>: output current during high state (source current)
- IoI: output current during low state (sink current)
- V<sub>i</sub>: allowed input voltage
- V<sub>ovp</sub>: overvoltage protection voltage
- V<sub>ih</sub>: input voltage that is interpreted as high level
- V<sub>il</sub>: input voltage that is interpreted as low level
- R<sub>i</sub>: input impedance value
- R<sub>o</sub>: output impedance value
- F<sub>co</sub>: input filter cut-off frequency
- INL: integral nonlinearity
- DNL: differential nonlinearity
- SINAD: signal-to-noise ratio
- Egain: gain error
- E<sub>offset</sub>: offset error in % full-scale-range (FSR)
- V<sub>offset</sub>: offset error
- V<sub>noise</sub>: DAC output noise
- $t_{settling}$ : output voltage settling time (transition time between 1/4 scale and 3/4 scale, including settling to  $\pm 0.024$  % within FSR)

# Sales partner / Contacts

# Headquarters

#### SmarAct GmbH

Schuette-Lanz-Strasse 9 26135 Oldenburg Germany

T: +49 441 – 800 87 90 Email: info-de@smaract.com www.smaract.com

## Japan

### **Physix Technology Inc.**

Ichikawa-Business-Plaza 4-2-5 Minami-yawata, Ichikawa-shi 272-0023 Chiba Japan T/F: +81 47 – 370 86 00

Email: info-jp@smaract.com www.physix-tech.com

### France

#### SmarAct GmbH

Schuette-Lanz-Strasse 9 26135 Oldenburg Germany

T: +49 441 – 80 08 79 956 Email: info-fr@smaract.com www.smaract.com

### **South Korea**

### **SEUM Tronics**

# 801, 1, Gasan digital 1-ro Geumcheon-gu Seoul, 08594, Korea

T: +82 2 868 – 10 02 Email: info-kr@smaract.com www.seumtronics.com

### Israel

#### Trico Israel Ltd.

P.O.Box 6172 46150 Herzeliya Israel

T: +972 9 – 950 60 74 Email: info-il@smaract.com www.trico.co.il

# **USA**

#### SmarAct Inc.

2140 Shattuck Ave. Suite 1103 Berkeley, CA 94704 United States of America

T: +1 415 – 766 9006 Email: info-us@smaract.com www.smaract.com